

Recovery Readiness Scoring System Using Wearable Data

Name: Yash Karsh

Date: 20/06/2025

1. Executive Summary

This report proposes a data-driven recovery readiness scoring system designed for integration with wearable health devices. It draws on scientific literature and real-world tools like WHOOP and Fitbit, incorporating metrics like HRV, RHR, sleep, and training load to guide post-exercise decision-making

2. Introduction

Post-Exercise recovery is a critical yet often overlooked component of physical performance and well-being. Without adequate recovery, individuals risk fatigue, injury, and reduced training effectiveness. In recent years, wearable technologies have empowered users to track physiological signals in real time – opening new opportunities for personalized fitness insights.

This report proposes a recovery readiness scoring system that leverages data from wearable devices to guide daily training decisions. Inspired tools such as WHOOP, Garmin Body Battery, and Fitbit Readiness Score, the proposed system integrates biometrics like Heart Rate Variability(HRV), Resting Heart Rate(RHR), sleep quality, and training intensity into a unified score. This score helps users determine whether to train hard, take it easy, or rest completely.

Backed by scientific research and grounded in practical deployment logic, the system is designed to be modular, evidence-based, and suitable for integration into digital wellness platforms.

3. Research & Assumptions

3.1 Existing Recovery Scoring Tools

Several fitness and wellness wearables use recovery scores to guide daily training readiness. Three leading tools are WHOOP, Garmin, and Fitbit:

WHOOP : Calculates a daily **Recovery Score (0–100%)** based on HRV, resting heart rate (RHR), sleep quality, respiratory rate, and skin temperature. It uses personalized baselines and presents scores in a color-coded format: Green (Ready), Yellow (Moderate), Red (Not Ready). WHOOP emphasizes recovery-driven training and has published internal studies showing fewer injuries and sustained performance when training is guided by their score.

Garmin: Offers a “Body Battery” metric derived from HRV and stress data. It tracks energy levels throughout the day and combines this with sleep and activity metrics. Garmin focuses on real-time energy monitoring rather than a single morning score.

Fitbit : The Daily Readiness Score (Premium feature) combines recent activity, HRV, and sleep. It classifies readiness into low, medium, or high and offers actionable suggestions like rest or training intensity

Metric	WHOOP	Garmin	Fitbit
HRV (Heart Rate Variability)	✓	✓	✓ (Premium)
RHR (Resting Heart Rate)	✓	✓	✓
Sleep Quality & Duration	✓	✓	✓
Respiratory Rate	✓	✓	✓
Skin Temperature	✓	□	✓ (Sense/Charge)
SpO ₂ (Blood Oxygen)	✓ (4.0 only)	✓	✓
Training Load / Strain	✓	✓ (Body Battery)	Limited
Stress Detection	□	✓	✓
Personalized Baseline Tracking	✓	✓	Limited

These platforms validate the need for a holistic, bio-signal-based recovery scoring methods.

3.2 Scientific Backing

The recovery score design is guided by the following scientifically validated physiological markers:

Heart Rate Variability (HRV): HRV measures the variation in time between heartbeats and reflects autonomic nervous system balance. A higher HRV generally indicates better recovery and resilience. Studies confirm HRV as a reliable indicator of training adaptation and overtraining risk (Plews et al., 2013; Buchheit, 2014).

Resting Heart Rate (RHR): Elevated RHR can indicate fatigue, illness, or incomplete recovery. It is often used in conjunction with HRV to assess training readiness.

Sleep Quality and Duration: Adequate sleep is essential for hormonal and muscular recovery. Poor sleep quality is linked to reduced HRV, higher RHR, and impaired performance (Fullagar et al., 2015).

Respiratory Rate & Skin Temperature: Used in WHOOP to detect early signs of illness or stress. Consistent increases can indicate reduced physiological resilience.

3.3 Assumptions

For the purposes of this scoring system, we make the following realistic assumptions:

User Demographics: Adults aged 25–45 with moderate to high fitness levels; gender-neutral.

Device Capabilities: The wearable can track HRV, RHR, sleep metrics, respiratory rate, and temperature.

Baseline Values: Each user's metrics are compared against their own rolling 30-day baseline to personalize the score.

Lifestyle Factors: Some metrics (like hydration, caffeine/alcohol use, and stress) may be collected via in-app self-report or passive tracking.

Note : "Currently available wearables like WHOOP 4.0, Fitbit Sense, and Oura Ring Gen 3 support the majority of required biometric signals. Where lifestyle data (e.g., hydration, stress) is not automatically sensed, user inputs via app-based journals can supplement the system."

4. Scoring System Design

4.1 Overview

The proposed recovery readiness system outputs a score on a 0–100 scale, which is color-coded for easy interpretation:

Green (70–100): Fully recovered – ready for intense training

Yellow (40–69): Partial recovery – light/moderate activity recommended

Red (0–39): Poor recovery – rest or very light activity advised

The score aggregates biometric and lifestyle data into a weighted average using a personalized baseline for each user.

4.2 Scoring Inputs

Below is the list of input metrics used in the score, including their purpose and suggested score range:

Metric	Description	Score Range	Data Source
Sleep Quality & Duration	Total hours, sleep stages, sleep efficiency	0–25	Wearable sleep tracker
HRV Deviation	Difference from personal baseline (higher = better)	0–25	Optical HRV sensor
Resting Heart Rate (RHR)	Elevation from baseline (lower = better)	0–15	Morning RHR value
Respiratory Rate	Nighttime rate vs. baseline (deviation indicates stress)	0–10	Wearable sensor
Skin Temperature	Deviation from baseline (e.g., +1°C = strain/illness)	0–10	Skin temp sensor
Training Load	Prior day's strain: duration × intensity	–15 to 0	App training log / sensor
Lifestyle Modifier	Hydration, stress, alcohol/caffeine impact	–15 to +5	User input or proxy

> Note: All metrics are normalized to a 0–100 scoring framework, then scaled using weightings described below.

4.3 Score Formula

A baseline sketch:

...

Recovery Score (0–100) =

$$\begin{aligned} &w1 \cdot \text{SleepScore} + w2 \cdot \text{HRV_deviation_score} + w3 \cdot \text{RHR_score} \\ &+ w4 \cdot \text{Resp_rate_score} + w5 \cdot \text{Temp_score} - w6 \cdot \text{TrainingLoad_penalty} \\ &+ w7 \cdot \text{Lifestyle_factor_adjustments} \end{aligned}$$

...

- Weightings (w1–w7) tuned via regression or domain-applied weights (e.g., sleep & HRV more heavily weighted).
- Baseline values and thresholds personalized (e.g. rolling 30-day average as baseline).:

$$\text{Recovery Score} = 0.25 \cdot \text{Sleep} + 0.25 \cdot \text{HRV} + 0.15 \cdot \text{RHR} + 0.10 \cdot \text{RespRate} + 0.10 \cdot \text{SkinTemp} + \text{TrainingLoadAdjustment} + \text{LifestyleModifier}$$

Where:

Sleep, HRV, RHR, etc. are individual sub-scores out of 100

Training Load Adjustment = –15 to 0 depending on strain (higher strain reduces score)

Lifestyle Modifier = –15 to +5 based on hydration, stress, alcohol, etc.

The weights reflect the physiological significance and scientific reliability of each signal. HRV and Sleep together account for 50% of the score due to their strong predictive power.

5. Example Personas

Persona 1: Alex — The Competitive Athlete

Age: 29

Fitness Level: High (training for a triathlon)

Sleep (last night): 8.2 hours, 90% efficiency

HRV: +10 ms above baseline

RHR: 3 bpm below baseline

Respiratory Rate: Normal (no deviation)

Skin Temperature: No deviation

Training Load (yesterday): High (intense cycling workout)

Lifestyle: Hydrated, no caffeine/alcohol, low stress

Scoring Breakdown:

Metric	Value	Sub-score (out of max)
Sleep	8.2 hrs, 90% eff	23/25
HRV Deviatioo	+10 ms	22/25
RHR	-3 bpm	13/15
Respiratory time	Normal	9/10
Skin Temperature	Normal	9/10
Training Load	High (-10)	-10
Lifestyle Identifier	+4	+4

Final Recovery Score:

$$= (0.25 \times 23) + (0.25 \times 22) + (0.15 \times 13) + (0.10 \times 9) + (0.10 \times 9) - 10 + 4$$
$$= 5.75 + 5.5 + 1.95 + 0.9 + 0.9 - 10 + 4$$
$$= 39.0 \rightarrow \text{Yellow Zone}$$

Interpretation: Alex is moderately recovered after a hard day. Light/moderate training is okay today.

Persona 2: Winnie — The Busy Beginner

Age: 35

Fitness Level: Low (office job, started jogging recently)

Sleep (last night): 5.5 hours, 75% efficiency

HRV: -20 ms below baseline

RHR: +6 bpm above baseline

Respiratory Rate: Elevated

Skin Temperature: +0.5°C above baseline

Training Load: Medium jog

Lifestyle: Stressed, 2 coffees, minimal water intake

Scoring Breakdown:

Metric	Value	Sub-score (out of max)
Sleep	5.5 hrs, 75% eff	13/25
HRV Deviatoo	-20 ms	10/25
RHR	+6 bpm	6/15
Respiratory time	Elevated	6/10
Skin Temperature	Slightly Elevated	6/10
Training Load	Medium (-5)	-5
Lifestyle Identifier	Poor (-10)	-10

Final Recovery Score:

$$= (0.25 \times 13) + (0.25 \times 10) + (0.15 \times 6) + (0.10 \times 6) + (0.10 \times 6) - 5 - 10$$

$$= 3.25 + 2.5 + 0.9 + 0.6 + 0.6 - 5 - 10$$

$$= -7.15 \rightarrow \text{Rounded to 0} \rightarrow \text{Red Zone}$$

Interpretation: Winnie is not recovered — high stress, low sleep, and HRV suppression signal a rest day is best.

6. System Implementation Proposal

6.1. System Overview

The recovery scoring system can be deployed within a mobile app or a wearable platform using the following architecture:

1. Input Layer

- Sensor Data: HRV, RHR, respiratory rate, skin temperature, sleep metrics from the wearable device
- User Inputs: Manual entries for hydration, stress levels, alcohol/caffeine intake, illness, etc.
- Training Log: Automatically detected or manually logged workouts (intensity \times duration)

2. Preprocessing Layer

- Baseline Normalization: Each user's daily inputs are compared against a personalized rolling 30-day average
- Noise Filtering: Remove outliers, smooth raw sensor data (especially HRV)
- Missing Data Handling: Use last-known-good values or estimate using trend smoothing

3. Score Computation Layer

* Apply the scoring formula:

$$\text{Recovery Score} = 0.25 \cdot \text{Sleep} + 0.25 \cdot \text{HRV} + 0.15 \cdot \text{RHR} + 0.10 \cdot \text{RespRate} + 0.10 \cdot \text{SkinTemp} + \text{Trai}$$


* Metrics are rescaled to standardized sub-scores before aggregation

4. Output Layer

- User Feedback:
- Recovery Score (0–100) with color coding
- Text-based advice (e.g., “Go hard today” or “Recovery day recommended”)
- Daily History View and Weekly Trends
- Optional push notifications (e.g., “Today’s score is low—prioritize rest”)

6.2. System Flow Diagram

Here’s a system block diagram of the proposed architecture:



> This architecture allows flexible integration with existing platforms like Garmin Connect, Fitbit Premium, or a custom fitness app.

✓ Optional Enhancements

- AI Layer: Use machine learning to optimize weightings per user
- Health Alerts: Notify users if abnormal trends (e.g., sustained elevated temperature) persist
- API Integration: Connect with other health data sources (Apple Health, Google Fit)

7. Conclusion

This report presented a comprehensive, data-driven recovery readiness scoring system designed for integration with wearable health technologies. By synthesizing scientific evidence and drawing inspiration from existing systems like WHOOP, Garmin, and Fitbit, the model leverages key physiological signals—such as heart rate variability, resting heart rate, sleep quality, and training load—to estimate recovery status on a 0–100 scale.

The system emphasizes personalization by using rolling baselines, and it accommodates lifestyle modifiers for stress, hydration, and other external factors. Practical implementation in a wearable app ecosystem was outlined, including input handling, signal preprocessing, and feedback delivery to users.

With modular design and evidence-backed logic, this system can guide users in making smarter training decisions, reducing the risk of overtraining or injury, and promoting sustainable fitness progress.

Future enhancements may include machine learning for adaptive weight tuning and third-party data integrations, making this scoring model both scalable and extensible in real-world digital health platforms.

8. Business Perspective

Market Potential of Health Monitoring Wearables

8.1 Industry Overview

The global health monitoring wearable devices market is experiencing rapid growth due to rising health awareness, increased adoption of fitness technologies, and the expansion of remote health monitoring. These devices go beyond step counting — they now offer insights into recovery, sleep, heart health, and stress, powered by AI and biosensors.

Market Size

According to Grand View Research, the global wearable healthcare market was valued at USD 26.8 billion in 2022, and it is expected to grow at a CAGR of 25.7% from 2023 to 2030.

Key Drivers

- Surge in lifestyle-related diseases
- Post-COVID interest in remote health diagnostics
- Expansion of health-conscious consumer base
- Growth in smart device ownership (phones + wearables)

8.2 Business Opportunities


a) B2C Market (Consumers)

- Fitness and wellness enthusiasts seeking personalized recovery/training insights
- Athletes and sportspersons looking to avoid injury and optimize performance
- Corporate wellness programs that incentivize employee health tracking

b) B2B Market (Businesses/Enterprises)

- Insurance Companies: Dynamic premium models based on wellness data
- Healthcare Providers: Integration into digital health records and remote patient monitoring (RPM)
- Gyms and Training Centers: Offer wearables with branded apps for client retention

8.3 Market Share of Major Wearable Device Brands



Apple and Fitbit dominate, especially in smartwatches and fitness bands, while other brands combined (including Xiaomi, Amazfit, Oura, boAt) hold a significant portion of the market.

- Emerging challengers: Xiaomi, Huawei, Garmin, Samsung, and niche players (e.g., Oura in smart rings).
- Market booming due to local brands (e.g., India's boAt) and niche devices (smart rings, medical bands).
- Overall market growing 4–10% CAGR, with smart rings a fast-growing niche

8.4 Monetization Models

- Subscription SaaS Model: Offer advanced insights, personalized recommendations, and coaching through a premium app
- Hardware + App Bundle: Sell the wearable with basic free app and optional premium upgrade
- White-label Licensing: License scoring system to existing health platforms or gyms

8.5 Risks & Challenges

- Privacy & Data Regulation (e.g., HIPAA, GDPR compliance)
- Sensor Accuracy & Validation
- User Retention: Maintaining long-term engagement with non-athlete

Conclusion

Recovery-based scoring systems are not only scientifically valid but commercially viable. As wearables continue shifting from activity tracking to *health insight engines*, businesses that invest in transparent, personalized, and engaging recovery analytics will gain a strong foothold in the \$100+ billion digital health market.


9. Bonus Insights : Simulated Data and Visual Analysis

To further illustrate the practicality and impact of the proposed recovery readiness system, this section includes simulated data visualizations, user interface mockups, and a system architecture diagram. These enhancements demonstrate how real-world users might interact with the scoring system and benefit from its insights.

10.1 Simulated Recovery Score Trend (5 Days)

A simulated user's recovery score over five days is shown in the chart below. The scores reflect daily variations based on different levels of sleep, HRV, and lifestyle inputs.

Figure: Simulated Recovery Score Over Time




This shows a clear dip mid-week followed by recovery on better-rested days.

10.2 Before vs. After Lifestyle Changes

The following graph compares recovery scores for a user before and after implementing healthier routines (improved hydration, consistent sleep, reduced caffeine). The improvement demonstrates the system's potential to promote behavior change.


Figure: Recovery Score Before vs After Lifestyle Changes



10.3 Visual Feedback Mockup

Below is a sample interface mockup showing how a user might receive their recovery score with color-coded guidance.

● *Mockup: In-App Recovery Score Display*



Message-driven feedback increases user engagement and supports daily decision-making.

10. References

Research Tools & Product Sources

- * WHOOP. (n.d.). *Understanding Recovery.* Retrieved from https://www.whoop.com
- * Fitbit. (2022). *What is the Daily Readiness Score?* Retrieved from [https://help.fitbit.com](https://help.fitbit.com/articles/en_US/Help_article/4203.htm)
- * Garmin. (n.d.). *Body Battery™ Energy Monitoring.* Retrieved from [https://www.garmin.com](https://www.garmin.com/en-US/blog/fitness/body-battery/)

Scientific Literature

- * Plews, D. J., Laursen, P. B., Stanley, J., Kilding, A. E., & Buchheit, M. (2013). *Training adaptation and heart rate variability in elite endurance athletes: Opening the door to effective monitoring.* Sports Medicine, 43(9), 773–781.
- * Buchheit, M. (2014). *Monitoring training status with HR measures: Do all roads lead to Rome?* Frontiers in Physiology, 5, 73.
- * Fullagar, H. H. K., et al. (2015). *Sleep and athletic performance: The effects of sleep loss on exercise performance, and physiological and cognitive responses to exercise.* Sports Medicine, 45(2), 161–186.

Optional White Papers & Validations

- * WHOOP Research Team. (2020). *The Impact of Recovery-Based Training on Performance and Injury Prevention*. Retrieved from https://www.whoop.com/the-locker